

The information in this document contains general descriptions of the technical options available, and may not always apply in individual cases. The required features should therefore be specified in each individual case at completion of contract.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information. Original images always lose a certain amount of detail when reproduced. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Global Business Unit Address

Siemens AG
Medical Solutions
Angiography, Fluoroscopic and
Radiographic Systems
Siemensstrasse 1
D-91301 Forchheim
Germany
Telephone: +49 9191 18-0
www.siemens.com/healthcare

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
D-80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
D-91052 Erlangen
Germany
Telephone: +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer

Siemens AG
Wittelsbacherplatz 2
D-80333 Muenchen
Germany

www.siemens.com/healthcare

MULTIX Swing

Fit for purpose. Easy to use.

www.siemens.com/healthcare

SIEMENS





SIEMENS

MULTIX Swing Pep up your routine.

MULTIX Swing lends pep to your X-ray room workflow and performance. As a compact, reliable, easy-to-use radiography solution with a broad applications spectrum, MULTIX Swing is a refreshing addition to the clinical environment.

MULTIX Swing minimizes training time, ensures smooth workflow and speeds up your patient throughput. The highly versatile system supports all orthopedic, thorax, as well as trauma and emergency diagnostic procedures – letting you perform your examinations with ease. MULTIX Swing is a cost-efficient, all-in-one X-ray room solution – one with the potential to bring new momentum into your imaging business.


German engineering for every budget

Easily accessible, the recessed foot switch enables efficient and safe operation.

The cassette is automatically centered upon insertion. This leads to reduced examination time.




All inclusive: floating tabletop, oscillating grid and dose measurement*. These features ensure ease of use, high image quality and conformance to radiation safety guidelines.



The floating tabletop with a weight capacity of up to 450 kg enables easy, comfortable examinations.



The synchronized tube and bucky tray movements ensure constant centering of the radiation beam on the cassette. This simplifies system positioning for faster examinations.




The tube can be rotated for cross-table exposures, expanding the spectrum of clinical applications.


* Option

Great user benefits, minimal space requirement

Chest radiography can be easily performed on a vertical chest stand.



The tube assembly enables vertical, oblique, horizontal and lateral acquisitions.



MULTIX Swing is specifically designed to provide a high degree of flexibility, even in a small room. The generator is integrated into the table and does not require additional space. This integration also facilitates and accelerates the installation.



In combination with a trolley, the wide applications range is further extended.



The unit is well suited to meet the needs of all patient types. Thanks to its ergonomic design and operator-friendly controls, patients can be examined comfortably.



Combined with the bucky wall stand VERTIX Swing (optional), the system offers you additional examination flexibility.

Life – our customer care solution

Increase your skills, technology and productivity

Life is the unique customer care solution from Siemens that helps you get the most from your investment.



Life. The continuous development of your ...

... skills:

On-site applications training. CD and Web-based e-learning tools to optimize your know-how to use the full potential of your system.

... productivity:

Siemens comprehensive portfolio of customized performance plans and variable shared service agreements for every type of service support as well as the Guardian Program™ for high system availability keep your equipment running at peak performance day in and day out.

... technology:

Keep your system state-of-the-art and expand your applications range with new functionalities for more versatility in use.

For more information about Life, go to:

www.siemens.com/life