



MammoDiagnost Mammography system



The **system**
you can **rely on**

PHILIPS



A smart solution for all applications

Introducing Philips MammoDiagnost, the highly versatile system for high volume screening, special views, magnification, and stereotactic biopsies. Its user-friendly and ergonomic design, perfect patient accessibility and superior image quality make the MammoDiagnost the system of choice for greater diagnostic confidence and bring the financial advantages associated with increased patient throughput.

And there's more...



Rise to the Challenge

In the field of medicine you face challenges every day. With a mammography system from Philips you can meet every challenge without even thinking about it. The MammoDiagnost is perfectly equipped for every application you might need. It is also extremely user and patient friendly and offers maximum comfort for both users and patients – there are numerous ways of positioning the system, and it has optimized compression and dose control.

Optimum dosage facility

Auto-selection of the best anode/filter combination using an optional unique four focal spot, dual-track X-ray tube ensures the lowest possible dose to suit specific breast characteristics – important for younger women with radiographically dense breasts. With W/Rh you can reduce the dosage by up to 60% compared to Mo/Mo exposures.



The possibility of upgrading to Digital Spot Imaging and digital stereotactics gives you flexibility when responding to future demands. The Philips' approach to digitization enables you to change over to digital radiography in stages and at a relatively low cost. In practice that means that investing in a MammoDiagnost will soon pay its way – through more efficient work processes, more effective results, and enhanced patient confidence.

Isocentric Rotation

Time-saving positioning of patients with motorized rotational movements, while maintaining the correct object height for all projections.



Pivoting Bucky's

Unique dual Bucky design (optional) for quick and easy switching between 18cm x 24cm and 24cm x 30cm formats; automatic pre-selection of collimation when object table is changed; manual replacement of either object table with stereotactic device or magnification table.



With its automatic compression optimization that adapts to individual breast characteristics and DoseWise dose control, MammoDiagnost is a system that takes patient safety and comfort to a whole new level.

Philips MammoDiagnost

- 25 years of experience in mammography
- The pioneer in digital mammography
- Technology you can trust: there are more than 5,000 Philips mammography units in medical facilities throughout the world


Philips – for the best in image quality

In order to be able to put your expertise to the best possible use for you and your patient you need images of outstanding quality. The smaller the finding, the greater the likelihood of a complete recovery. With Philips you have a partner, who offers the best in image quality – to help you build trust with your patients and for an easier and more accurate diagnosis. Because timely detection is everything.



Philips MammoDiagnost

- Outstanding image quality, from the pioneer of digital X-ray
- Sophisticated design for ideal ergonomics
- CUSTOMerCare life cycle support packages
- Modular system family for customized setups
- Optimized initial lifetime costs



**Philips Medical Systems is part
of Royal Philips Electronics**

Interested?

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be glad to hear from you.

On the web

www.medical.philips.com

Via e-mail

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems

Global Information Center

P.O.Box 1168

5602 BD Eindhoven

The Netherlands

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +31 40 27 63005

Latin America

Tel: +55 11 2125 0764

North America

Tel: +1 800 229 6417

© Koninklijke Philips Electronics N.V. 2004

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems DMC GmbH reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.

4522 981 91891/714* AUG 2004

